



Peninsula Laboratories, LLC
A Member of the Bachem Group
305 Old County Road, San Carlos, CA 94070
Tel: (800) 922-1516 • (650) 592-5392
Fax: (650) 595-4071
www.bachem.com

Cat. No. S-4020
Dynorphin A (1-8)

Specifications

Number of Slides:	150
Chromogens:	DAB (brown), Ni-intensified DAB (dark brown), and AEC (red)
Sample:	Paraffin-embedded or frozen sections
Shelf Life:	One year from production date
Storage:	4°C
Assay Time:	2-3 hours following primary incubation

***CROSS REACTION WITH:**

Dynorphin A (1-8)
Prodynorphin (209-240) (porcine)
Dynorphin A
Dynorphin A (1-13)
Dynorphin B
Met-Enkephalin
β-Endorphin (porcine)
α-Neoendorphin

%:

100
1
2
2
0
0
0
0

Kit Components

- 30% Hydrogen peroxide Solution
- Blocking Solution
- Primary Antibody
- Biotinylated Secondary Antibody
- SA-HRP Conjugate Solution
- DAB Concentrated Solution
- Nickel-intensified DAB Concentrated Solution
- AEC Concentrated Solution
- Substrate Concentrated Buffer
- Hydrogen Peroxide Concentrate
- Hematoxylin Counterstain
- Mounting Solution

***tested by RIA**