



Peninsula Laboratories, LLC
A Member of the Bachem Group
305 Old County Road, San Carlos, CA 94070
Tel: (800) 922-1516 • (650) 592-5392
Fax: (650) 595-4071
www.bachem.com

Cat. No. S-4024
Neurokinin A

Specifications

Number of Slides: 150
Chromogens: DAB (brown), Ni-intensified DAB (dark brown), and AEC (red)
Sample: Paraffin-embedded or frozen sections
Shelf Life: One year from production date
Storage: 4°C
Assay Time: 2-3 hours following primary incubation

***CROSS REACTION WITH:**

Neurokinin A
Kassinin
Neurokinin K (porcine)
Neurokinin B
Substance P
Physalaemin
Neuropeptide Y (porcine)
Neuromedin B (porcine)
ACTH (1-39) (human)
[Arg⁸]-Vasopressin

%:

100
100
100
80
0
0
0
0
0
0
0

Kit Components

- 30% Hydrogen peroxide Solution
- Blocking Solution
- Primary Antibody
- Biotinylated Secondary Antibody
- SA-HRP Conjugate Solution
- DAB Concentrated Solution
- Nickel-intensified DAB Concentrated Solution
- AEC Concentrated Solution
- Substrate Concentrated Buffer
- Hydrogen Peroxide Concentrate
- Hematoxylin Counterstain
- Mounting Solution

***tested by RIA**