

# Peninsula Laboratories, LLC

A Member of the Bachem Group 305 Old County Road, San Carlos, CA 94070 Tel: (800) 922-1516 • (650) 592-5392 Fax: (650) 595-4071 www.bachem.com

# Cat. No. S-4034 Kisspeptin-13 (4-13) (human)

## **Specifications**

| Number of Slides:<br>Chromogens: | 150<br>DAB (brown), Ni-intensified DAB (dark brown), and AEC<br>(red) |
|----------------------------------|---|
| Sample:                          | Paraffin-embedded or frozen sections                                  |
| Shelf Life:                      | One year from production date   |
| Storage:                         | 4°C   |
| Assay Time:                      | 2-3 hours following primary incubation                                |

%:

100

0

0

# \*CROSS REACTION WITH:

| Kisspeptin-13 (4-13) (human)  |  |
|-------------------------------|--|
| Kisspeptin-54 (human)         |  |
| Kisspeptin-54 (27-54) (human) |  |
| Antho-Rwamide I               |  |
| Antho-RWamide II              |  |
|                               |  |

#### \*tested by RIA

## Kit Components

- •30% Hydrogen peroxide Solution
- •Blocking Solution
- •Primary Antibody
  - •Biotinylated Secondary Antibody
  - •SA-HRP Conjugate Solution
  - •DAB Concentrated Solution
    - •Nickel-intensified DAB Concentrated Solution
    - •AEC Concentrated Solution
    - •Substrate Concentrated Buffer
    - •Hydrogen Peroxide Concentrate
    - •Hematoxylin Counterstain
    - Mounting Solution