

Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070 Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071 www.bachem.com

FITC Labeled Monoclonal Antibody To Mouse MHC Class II

Monoclonal antibody ER-TR 3 is one member of a family of monoclonal antibodies (ER-TR 3, ER-TR 2, ER-TR 1) which detect MHC class II antigens encoded by the murine Ia region of the H-2 complex, corresponding to the human HLA-DR region. They are valuable tools for studying T helper cell interaction with class II positive antigen presenting cells (dendritic cells, B-cells, macrophages). These antibodies also offer new possibilities for studying the development of T helper cells since they also stain stromal cells in the thymus.

Product Number:T-2117Clone:ER-TR3Host species, isotype:Rat IgG2bQuantity:100μg

Format: Affinity purified, FITC labelled, liquid

Supplied as 1ml solution. This stock solution contains 0.1mg/ml IgG, phosphate buffered saline pH 7.2 (PBS), 10mg/ml bovine serum albumin (BSA) as a stabilizer and 0.09% sodium azide

as a preservative.

Stability: Original vial: 6 months at 4° - 8°C

Applications: Has been described to work in FACS.

Approximate working dilution:

Optimal dilutions should be determined by the end user.

Suggested positive control: mouse spleen.

Immunogen: Murine thymic reticulum.

Antigen, epitope: MHC Class II antigens are heterodimers consisting of one α -

chain(31-34kD) and one β -chain (26-29kD).

Antigen distribution Isolated cellsThe antigen is found on dendritic cells, B-cells

and macrophages. The level of antigen detected by ER-TR1, ER-TR 2 and ER-TR 3 differs from strain to strain (see table

below).

Tissue Sections: The antigen is found on B-cells,

interdigitating cells and macrophages in peripheral lymphoid organs but is absent from T-cells. It is also expressed as a fine reticular pattern on stromal thymic cells of the cortex and as a confluent pattern on stromal thymic cells of the medulla.



Peninsula Laboratories, LLC

A Member of the Bachem Group

305 Old County Road, San Carlos, CA 94070 Tel: (800) 922-1516 • (650) 592-5392

Fax: (650) 595-4071 www.bachem.com

Distribution of ER-TR 1, ER-TR 2 and ER-TR 3 among mouse strains with independent and recombinant haplotypes*

Strain	Haplotype							Clone		
	K	Α	В	J	Ε	С	D	ER-TR1	ER-TR2	ER-TR3
C3H/HeJ	k	k	k	k	k	k	k	48*	46	46
AKR	k	k	k	k	k	k	k	54	52	54
B10.BR	k	k	k	k	k	k	k	59	58	62
B10.ScSn	b	b	b	b	b	b	b	4	5	50
Balb/b	b	b	b	b	b	b	b	4	3	39
B10.D2/n	d	d	d	d	d	d	d	56	5	54
Balb/c	d	d	d	d	d	d	d	45	3	44
DBA/2	d	d	d	d	d	d	d	27	4	47
B10.G	q	q	q	q	q	q	q	53	4	46
DBA/1	q	q	q	q	q	q	q	52	6	54
SWR/J	q	q	q	q	q	q	q	49	3	49
A.SW	S	S	S	S	S	S	S	4	20	6
B10.M	f	f	f	f	f	f	f	4	5	3
B10.RIII	r	r	r	r	r	r	r	39	39	40
B10.AQR	q	k	k	k	k	d	d	52	52	51
B10.T(6R)	q	q	q	q	q	q	d	50	3	52
A.TL	S	k	k	k	k	k	d	29	52	51
A.TH	S	S	S	S	S	S	d	5	49	7

^{*} Percentage of labelled cells, determined by FACS analysis of spleen cell suspensions

Specificity: Mouse: cells expressing MHC class II antigens

Other: negative on human, other unknown

Selected references

Van Vliet, E., et al.: Monoclonal Antibodies to Stromal Cell Types of the Mouse Thymus. Eur. J. Immunol. <u>14</u>, 524-529 (1984)

Van Vliet, E., et al.: Stromal Cell Types in the Developing Thymus of the Normal and Nude Mouse Embryo. Eur. J. Immunol. <u>15</u>, 675-681 (1985)

For in vitro research only. Caution: this product contains sodium azide, a poisonous and hazardous substance.