



PENINSULA LABORATORIES, LLC

A MEMBER OF THE BACHEM GROUP

305 OLD COUNTY ROAD, SAN CARLOS, CALIFORNIA 94070

(800) 922-1516 * (650) 592-5392 * FAX (650) 595 4071

www.bachem.com

T-4005

Rabbit anti-Angiotensin II IgG

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
Angiotensin II	100
[Val ⁵]-Angiotensin II	100
Angiotensin III	100
Renin Substrate	0.9
Angiotensin I	0.5
ANF (1-28) (human)	0
[Arg ⁸]-Vasopressin	0
[Sar ¹ ,Ile ⁸] Angiotensin II	0