

## PENINSULA LABORATORIES, LLC

A Member of the Bachem Group 305 Old County Road, San Carlos, CA 94070 Tel: (800) 922-1516 • (650) 592-5392 Fax: (650) 595-4071 www.bachem.com

## Cat. No. **T-4246** Rabbit anti-VIP (human, bovine, porcine, rat)

## **Specifications:**

Preparation:	This polyclonal antibody was collected from rabbits immunized with a synthetic peptide as the immunogen.	
Amount:	Each vial contains 50 ul of whole serum in lyophilized form. Reconstitute with 50 ul of distilled water for the equivalent of undiluted antiserum.	
Storage:	This antibody is supplied as a lyophilized powder. Store at 2° to 8°C upon receipt. For long-term storage, store lyophilized powder at -20°C. Once reconstituted, the materials should be stored frozen in aliquot at -20°C to -70°C. <b>Avoid repeated</b> <b>freeze/thaw cycles.</b>	
Applications:	Immunohistochemistry, Immunocytochemistry, and other applications.	
Sample Types:	The antibody can react with a variety of sample preparation including	
	frozen tissue sections, paraffin embedded tissue sections and fixed cell cultures.	
Recommended starting dilution for Immunohistochemical staining:		

Immunofluorescence: 1:250 to 1:16000 PAP or ABC: 1:250 to 1:16000 LAB-SA: 1:250 to 1:16000 It is advised for each laboratory to determine optimal dilutions for a specific application.

## **Cross Reaction:**

Peptide	%
VIP (human, porcine, rat)	100
VIP (guinea pig)	100
VIP 10-28 (human, porcine, rat)	100
VIP (chicken)	<]
PACAP-27 (human, ovine, rat)	< 0.02
VIP 1-12 (human, porcine, rat)	0
PHM-27 (human)	0
Substance P	0
Endothelin-1 (human, porcine, canine, rat, mouse, bovine)	0
Secretin (porcine)	0
Glucagon (human)	0
Galanin (rat)	0
Somatostatin	

Cross-Reactivity data is taken from Radioimmunoassay applications