

Catalog # CEK-011-10

Standard Size 96 Wells

Sensitivity Minimum Detectable Concentration = 12 pg/ml

Range 0 - 10000 pg/ml

Linear Range 12 - 535 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	BNP-32 (Porcine)	100
	BNP-26 (Porcine)	100
	BNP-32 (Human)	0
	BNP-32 (Rat)	0
	BNP-45 (Rat)	0
	Alpha-ANP (1-28) (Human, Canine, Ovine)	0

Sequence Ser - Pro - Lys - Thr - Met - Arg - Asp - Ser - Gly - Cys - Phe - Gly - Arg - Arg - Leu - Asp - Arg - Ile - Gly - Ser - Leu - Ser - Gly - Leu - Gly - Cys - Asn - Val - Leu - Arg - Arg - Tyr

Disulfide Bridge: Cys10 - Cys26

Standard Curve Image

