

Catalog # EK-011-10

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.45 ng/ml

Range 0-100 ng/ml

Linear Range 0.45-6.3 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	BNP-32 (Porcine)	100
	BNP-26 (Porcine)	100
	BNP-32 (Human)	0
	BNP-32 (Rat)	0
	BNP-45 (Rat)	0
	α-ANP (1-28) (Human, Canine)	0

Sequence Ser - Pro - Lys - Thr - Met - Arg - Asp - Ser - Gly - Cys - Phe - Gly - Arg - Arg - Leu - Asp - Arg - Ile - Gly - Ser - Leu - Ser - Gly - Leu - Gly - Cys - Asn - Val - Leu - Arg - Arg - Tyr

Disulfide Bridge: Cys10 - Cys26

Image

