

Catalog # EK-022-06

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.15 ng/ml

Range 0-100 ng/ml

Linear Range 0.15-2.23 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	β-Endorphin (Camel, Bovine, Ovine)	100
	β-Endorphin (Porcine)	100
	β-Endorphin (Rat)	100
	β-Endorphin (Equine)	100
	α-Endorphin	100
	β-Endorphin (Human)	40
	Met-Enkephalin	0
	Leu-Enkephalin	0
	α-MSH	0
	ACTH (Human)	0
	PACAP-38 (Human, Ovine, Rat)	0

Sequence Tyr - Gly - Gly - Phe - Met - Thr - Ser - Glu - Lys - Ser - Gln - Thr - Pro - Leu - Val - Thr - Leu - Phe - Lys - Asn - Ala - Ile - Ile - Lys - Asn - Ala - His - Lys - Lys - Gly - Gln

Image

