

Catalog # EK-022-32

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.13ng/ml

Range 0-100ng/ml

Linear Range 0.13-4.46 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	β-Endorphin (Porcine)	100
	β-Endorphin (Rat)	100
	β-Endorphin (Camel, Ovine, Bovine)	100
	β-Endorphin (Equine)	100
	α-Endorphin	100
	r-Endorphin	60
	β-Endorphin (Human)	50
	Met-Enkephalin	0
	Leu-Enkephalin	0

Sequence Tyr - Gly - Gly - Phe - Met - Thr - Ser - Glu - Lys - Ser - Gln
 - Thr - Pro - Leu - Val - Thr - Leu - Phe - Lys - Asn - Ala - Ile
 - Val - Lys - Asn - Ala - His - Lys - Lys - Gly - Gln

Standard Curve Image

