

Catalog #	EK-060-04														
Standard Size	96 wells														
Sensitivity	Minimum Detectable Concentration = 0.31 ng/ml														
Range	0-100 ng/ml														
Linear Range	0.31 – 2.41ng/ml														
Cross Reactivity	<table> <tr> <th>Peptide</th><th>%Cross-reactivity</th></tr> <tr> <td>Cortistatin-14 (Rat)</td><td>100</td></tr> <tr> <td>Cortistatin-29 (Rat)</td><td>100</td></tr> <tr> <td>Cortistatin-29 (Human)</td><td>100</td></tr> <tr> <td>Cortistatin-17 (Human)</td><td>100</td></tr> <tr> <td>Somatostatin</td><td>100</td></tr> <tr> <td>Somatostatin-28</td><td>100</td></tr> </table>	Peptide	%Cross-reactivity	Cortistatin-14 (Rat)	100	Cortistatin-29 (Rat)	100	Cortistatin-29 (Human)	100	Cortistatin-17 (Human)	100	Somatostatin	100	Somatostatin-28	100
Peptide	%Cross-reactivity														
Cortistatin-14 (Rat)	100														
Cortistatin-29 (Rat)	100														
Cortistatin-29 (Human)	100														
Cortistatin-17 (Human)	100														
Somatostatin	100														
Somatostatin-28	100														
Sequence	Pro - Cys - Lys - Asn - Phe - Phe - Trp - Lys - Thr - Phe - Ser - Ser - Cys - Lys														

Standard Curve Image

