

Catalog #	EK-060-12										
Standard Size	96 wells										
Sensitivity	Minimum Detectable Concentration = 0.1ng/ml										
Range	0-100ng/ml										
Linear Range	0.1 – 0.97ng/ml										
Cross Reactivity	<table> <tr> <th>Peptide</th><th>%Cross-reactivity</th></tr> <tr> <td>Cortistatin-29 (Human)</td><td>100</td></tr> <tr> <td>Cortistatin-29 (Rat)</td><td>100</td></tr> <tr> <td>Cortistatin-14 (Rat)</td><td>100</td></tr> <tr> <td>Cortistatin-17 (Human)</td><td>100</td></tr> </table>	Peptide	%Cross-reactivity	Cortistatin-29 (Human)	100	Cortistatin-29 (Rat)	100	Cortistatin-14 (Rat)	100	Cortistatin-17 (Human)	100
Peptide	%Cross-reactivity										
Cortistatin-29 (Human)	100										
Cortistatin-29 (Rat)	100										
Cortistatin-14 (Rat)	100										
Cortistatin-17 (Human)	100										

Sequence pGlu - Glu - Gly - Ala - Pro - Pro - Gln - Gln - Ser - Ala - Arg - Arg - Asp - Arg - Met - Pro - Cys - Arg - Asn - Phe - Phe - Trp - Lys - Thr - Phe - Ser - Ser - Cys - Lys

Disulfide Bridge: Cys17 - Cys28

Standard Curve Image

