

Catalog # EK-071-05

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.06 ng/ml

Range 0-100 ng/ml

Linear Range 0.06-1.48 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Urotensin II (Human)	100
	Urotensin II(5-11) (Human)	22
	Urotensin II (Rat)	18.2
	Urotensin II Related Peptide (Human,Mouse, Rat)	15.7
	Urotensin II (Mouse)	15.1
	Urotentin II (4-11) (Human)	15
	Urotensin II (102-125)-Prepro(Human)	0
	Urotensin I (Catostomus commersoni)	0
	Endothelin-1	0
	Angiotensin II	0
	ADM (Rat)	0
	PAMP-20 (Human)	0
	α -ANP-28 (Human)	0
	CGRP (Rat)	0
	CNP-22	0
	BNP-25 (Rat)	0
	BNP-32 (Rat)	0
	Bradykinin	0
	Neurotensin	0

Sequence Glu - Thr - Pro - Asp - Cys - Phe - Trp - Lys - Tyr - Cys - Val

Disulfide Bridge: Cys5 - Cys10

Standard Curve
Image

