

Catalog #	RK-002-01	
Standard Size	125 RIA tubes	
Radioactivity	This kit contains 1.5 µCi of Iodine 125	
Range	1-128 pg/tube of 100 µl Standard Solution = 10-1280 pg/ml	
Lowest Detection Limit	36.6 pg/ml	
Cross Reactivity	<i>Peptide</i>	<i>% Crossreactivity</i>
	Angiotensin I (Human)	100
	[Asn ¹ , Val ⁵ , Asn ⁹]-Angiotensin I (Salmon)	100
	[Val ¹ , Asn ⁹]-Angiotensin-I (Bullfrog)	100
	[Val ¹ , Ser ⁹]-Angiotensin-I (Fowl)	100
	[Asn ¹ , Val ⁵]-Angiotensin-I (Goosefish)	13
	Angiotensin II (Human)	0
	Angiotensin III (Human)	0
	Neurotensin	0
	Substance P	0
	[Arg ⁸]-Vasopressin	0
	Endothelin-1 (Human, Canine, Porcine, Rat, Mouse)	0
	BNP-32 (Human)	0
	α-ANP (1-28) (Human, Canine)	0
	Renin Substrate (Human)	0
Sequence	Asp - Arg - Val - Tyr - Ile - His - Pro - Phe - His - Leu	