

Catalog #	RK-003-01
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	50-6400 pg/tube of 100µl Standard Solution = 500-64000 pg/ml
Lowest Detection Limit	2174 pg/ml
Cross Reactivity	not canine
Sequence	Ser - Cys - His - Leu - Pro - Trp - Ala - Ser - Gly - Leu - Glu - Thr - Leu - Asp - Ser - Leu - Gly - Gly - Val - Leu - Glu - Ala - Ser - Gly - Tyr - Ser - Thr - Glu - Val - Val - Ala - Leu - Ser - Arg - Leu - Gln - Gly - Ser - Leu - Gln - Asp - Met - Leu - Trp - Gln - Leu - Asp - Leu - Ser - Pro - Gly - Cys