

Catalog # RK-005-14

Standard Size 125 RIA tubes

Radioactivity This kit contains 1.5 µCi of Iodine 125

Range 1-128pg/tube in 100 µl of Standard Solution = 10-1280pg/ml

Lowest Detection Limit 67.4 pg/ml

Cross Reactivity	<i>Peptide</i>	<i>% Cross-reactivity</i>
	β-ANP (Human)	100
	α-ANP (1-28) (Human, Canine)	3
	Urodilatin	3
	ANP (Rat)	0
	BNP-32 (Human)	0
	CNP-22 (Porcine)	0

Sequence Ser - Leu - Arg - Arg - Ser - Ser - Cys - Phe - Gly - Gly - Arg - Met -
Asp - Arg - Ile - Gly - Ala - Gln - Ser - Gly - Leu - Gly - Cys - Asn -
Ser - Phe - Arg - Tyr

Disulfide Bridge: Cys7 - Cys51 and Cys23 - Cys35