

Catalog #	RK-005-18
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube of 100ul Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	38.7 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Urodilatin	100
	α-ANP (1-28) (Human, Canine)	100
	ANP (Rat)	100
	β-ANP (Human)	100
	BNP-32 (Human)	0
	CNP-22 (Porcine)	0
	Endothelin-1 (Human, Porcine, Rat, Mouse, Canine, Bovine)	0

Sequence Thr - Ala - Pro - Arg - Ser - Leu - Arg - Arg - Ser - Ser - Cys - Phe - Gly - Gly - Arg - Met - Asp - Arg - Ile - Gly - Ala - Gln - Ser - Gly - Leu - Gly - Cys - Asn - Ser - Phe - Arg - Tyr

Disulfide Bridge: Cys11 - Cys27