

Catalog #	RK-005-19
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube in 100 µl of Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	39 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Prepro-ANP (26-55) (Human)	100
	α-ANP (1-28) (Human, Canine, Porcine, Bovine)	0
	ANP (Rat)	0
	β-ANP (Human)	0
	Prepro-ANP (56-92) (Human)	0
	BNP-32 (Human)	0
	CNP-22 (Human, Porcine, Rat)	0
	Insulin	0
	IGF-1	0

Sequence Asn - Pro - Met - Tyr - Asn - Ala - Val - Ser - Asn - Ala - Asp - Leu - Met - Asp - Phe - Lys - Asn - Leu - Leu - Asp - His - Leu - Glu - Glu - Lys - Met - Pro - Leu - Glu - Asp