

Catalog # RK-005-20

Standard Size 125 RIA tubes

Radioactivity This kit contains 1.5 µCi of Iodine 125

Range 10-1280 pg/tube of 100 µl Standard Solution or
100-12,800 pg/ml

**Lowest
Detection Limit** 175.3 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Prepro-ANP (56-92) (Human)	100
	Prepro-ANP (26-92) (Human)	100
	α-ANP (1-28) (Human, Canine)	0
	β-ANP (Human)	0
	[Arg ⁸]-Vasopressin	0
	Oxytocin	0
	Somatostatin	0

Sequence Glu - Val - Val - Pro - Pro - Gln - Val - Leu - Ser - Glu - Pro
- Asn - Glu - Glu - Ala - Gly - Ala - Ala - Leu - Ser - Pro - Leu
- Pro - Glu - Val - Pro - Pro - Trp - Thr - Gly - Glu - Val - Ser
- Pro - Ala - Gln - Arg