

Catalog #	RK-010-08
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube in 100 µl of Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	45.2 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Adrenomedullin (Rat)	100
	Adrenomedullin (Human)	0
	Pro-Adrenomedullin N20 (Human)	0
	Pro-Adrenomedullin N20 (Rat)	0
	Amylin (Human)	0
	ET-1 (Human, Bovine, Canine, Mouse, Porcine, Rat)	0
	CGRP (Rat)	0
	Calcitonin (Rat)	0
	Calcitonin (Salmon)	0

Sequence Tyr - Arg - Gln - Ser - Met - Asn - Gln - Gly - Ser - Arg - Ser - Thr - Gly - Cys - Arg - Phe - Gly - Thr - Cys - Thr - Met - Gln - Lys - Leu - Ala - His - Gln - Ile - Tyr - Gln - Phe - Thr - Asp - Lys - Asp - Lys - Asp - Gly - Met - Ala - Pro - Arg - Asn - Lys - Ile - Ser - Pro - Gln - Gly - Tyr - NH₂

Disulfide Bridge: Cys14 - Cys19