

Catalog #	RK-010-10				
Standard Size	125 RIA tubes				
Radioactivity	This kit contains 1.5 µCi of Iodine 125				
Range	1-128 pg/tube in 100 µl of Standard Solution = 10~1280 pg/ml				
Lowest Detection Limit	49.81 pg/ml				
Cross Reactivity	<table> <tr> <th><i>Peptide</i></th><th><i>%Cross-reactivity</i></th></tr> <tr> <td>Pro-Adrenomedullin N-20 (Rat)</td><td>100</td></tr> </table>	<i>Peptide</i>	<i>%Cross-reactivity</i>	Pro-Adrenomedullin N-20 (Rat)	100
<i>Peptide</i>	<i>%Cross-reactivity</i>				
Pro-Adrenomedullin N-20 (Rat)	100				
Sequence	Ala - Arg - Leu - Asp - Thr - Ser - Ser - Gln - Phe - Arg - Lys - Lys - Trp - Asn - Lys - Trp - Ala - Leu - Ser - Arg - NH2				