

Catalog #	RK-010-52
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	2-256 pg/tube in 100µl of Standard Solution = 20-2560 pg/ml
Lowest Detection Limit	76 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Intermedin/Adrenomedullin-2 (Rat)	100
	Intermedin/Adrenomedullin-2 (Human)	100
	Adrenomedullin (Rat)	0
	Adrenomedullin (Human)	0
	CGRP (Rat)	0
	CGRP (Human)	0

Sequence Pro - His - Ala - Gln - Leu - Leu - Arg - Val - Gly - Cys - Val - Leu - Gly - Thr - Cys - Gln - Val - Gln - Asn - Leu - Ser - His - Arg - Leu - Trp - Gln - Leu - Val - Arg - Pro - Ser - Gly - Arg - Arg - Asp - Ser - Ala - Pro - Val - Asp - Pro - Ser - Ser - Pro - His - Ser - Tyr - NH₂ - Disulfide bridge between Cys10 & Cys15