

Catalog #	RK-010-62
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	20-2560 pg/tube of 100µl standard solution = 200-25600 pg/ml
Lowest Detection Limit	246 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
-------------------------	----------------	---------------------------

Intermedin /AM-2, Prepro-(57-92) (Human)	100
Intermedin / AM-2, Prepro-(25-66) (Human)	0
Intermedin / AM-2 (Human)	0
Intermedin-53 (Human)	0
Intermedin / AM-2 (8-47) (Human)	0
Adrenomedullin / AM-1 / ADM-1 (Human)	0
Adrenomedullin / AM-1 / ADM-1 (Rat)	0
Prepro-Adrenomedullin N-20 (Human) (PAMP-20)	0
Adrenomedullin, Pro-(45-92) (Human)	0
Adrenomedullin, Pro-(153-185) (Human)	0
CGRP (Human)	0

Sequence	Pro - Val - Val - Trp - Lys - Leu - His - Arg - Ala - Leu - Gln - Ala - Gln - Arg - Gly - Ala - Gly - Leu - Ala - Pro - Val - Met - Gly - Gln - Pro - Leu - Arg - Asp - Gly - Gly - Arg - Gln - His - Ser - Gly - Pro
-----------------	---