

Catalog # RK-011-03

Standard Size 125 RIA tubes

Radioactivity This kit contains 1.5 µCi of Iodine 125

Range 1-128 pg/tube of 100µl Standard Solution = 10-1280 pg/ml

Lowest Detection Limit 36 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	BNP-32 (Human)	100
	BNP-32 (Rat)	0
	BNP-45 (Rat)	0
	BNP-26 (Porcine)	0
	Prepro BNP (1-46) (Human)	0
	Prepro BNP (22-46) (Human)	0
	α-ANP (1-28) (Human, Canine)	0
	Endothelin-1 (Human, Rat, Porcine, Mouse, Canine, Bovine)	0
	Angiotensin II (Human)	0
	VIP (Human, Rat, Porcine)	0

Sequence Ser - Pro - Lys - Met - Val - Gln - Gly - Ser - Gly - Cys - Phe - Gly - Arg - Lys - Met - Asp - Arg - Ile - Ser - Ser - Ser - Ser - Gly - Leu - Gly - Cys - Lys - Val - Leu - Arg - Arg - His – Disulfide Bridge between Cys10 & Cys26

Standard Curve Image

