

Catalog #	RK-011-06	
Standard Size	125 RIA tubes	
Radioactivity	This kit contains 1.5 µCi of Iodine 125	
Range	1-128 pg/tube in 100µl of Standard Solution = 10-1280 pg/ml	
Lowest Detection Limit	65.52 pg/ml	
Cross Reactivity	Peptide	%Cross-reactivity
	Pro-BNP (22-46) (Human)	100
	BNP-32 (Human)	0
	BNP-32 (Rat)	0
	BNP-32 (Porcine)	0
	BNP-26 (Porcine)	0
	α-ANP (1-28) (Human, Canine)	0
	Endothelin-1 (Human, Rat, Porcine, Mouse, Canine, Bovine)	0
	CNP-22 (Human, Porcine, Rat)	0
Sequence	Asn - His - Leu - Gln - Gly - Lys - Leu - Ser - Glu - Leu - Gln - Val - Glu - Gln - Thr - Ser - Leu - Glu - Pro - Leu - Gln - Glu - Ser - Pro - Arg	

Standard Curve Image

