

Catalog #	RK-011-08
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube in 100µl of Standard Solution =10-1280 pg/ml
Lowest Detection Limit	50.85 pg/ml

Cross Reactivity	<i>Peptide</i>	<i>%Cross-reactivity</i>
	BNP-26 (Porcine)	100
	BNP-32 (Porcine)	48
	BNP-32 (Canine)	42
	BNP-32 (Human)	0
	BNP-32 (Rat)	0
	α-ANP (1-28) (Human, Canine)	0
	β-ANP (Human)	0
	ANP (Rat)	0
	Substance P	0

Sequence	Asp - Ser - Gly - Cys - Phe - Gly - Arg - Arg - Leu - Asp - Arg - Ile - Gly - Ser - Leu - Ser - Gly - Leu - Gly - Cys - Asn - Val - Leu - Arg - Arg - Tyr ~ Disulfide bridge between Cys4 & Cys20
-----------------	---