

Catalog #	RK-011-14	
Standard Size	125 RIA tubes	
Radioactivity	This kit contains 1.5 µCi of Iodine 125	
Range	1-128 pg/tube in 100µl of Standard Solution = 10-1280 pg/ml	
Lowest Detection Limit	56.7 pg/ml	
Cross Reactivity	<i>Peptide:</i>	<i>%Cross-reactivity</i>
	BNP-32 (Rat)	100
	BNP-45 (Rat)	100
	BNP-26 (Porcine)	0
	ANP (Rat)	0
	β-ANP (Human)	0
	Endothelin-1	
	(Human, Rat, Porcine, Mouse, Canine, Bovine)	0
	α-ANP (1-28) (Human, Canine)	0
Sequence	Asn - Ser - Lys - Met - Ala - His - Ser - Ser - Ser - Cys - Phe - Gly - Gln - Lys - Ile - Asp - Arg - Ile - Gly - Ala - Val - Ser - Arg - Leu - Gly - Cys - Asp - Gly - Leu - Arg - Leu - Phe ~ Disulfide bridge between Cys10 & Cys26	