

Catalog #	RK-011-24
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	10-1280 pg/tube in 100µl of Standard Solution = 100-12,800 pg/ml
Lowest Detection Limit	250 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Pro-BNP (1-46) (Human)	100

Sequence

His - Pro - Leu - Gly - Ser - Pro - Gly - Ser - Ala - Ser - Asp - Leu - Glu - Thr - Ser - Gly - Leu - Gln - Glu - Gln - Arg - Asn - His - Leu - Gln - Gly - Lys - Leu - Ser - Glu - Leu - Gln - Val - Glu - Gln - Thr - Ser - Leu - Glu - Pro - Leu - Gln - Glu - Ser - Pro - Arg

Standard Curve Image

