

Catalog #	RK-011-26	
Standard Size	125 RIA tubes	
Radioactivity	This kit contains 1.5 µCi of Iodine 125	
Range	10-1280 pg/tube in 100µl of Standard Solution = 100-12,800 pg/ml	
Lowest Detection Limit	144 pg/ml	
Cross Reactivity	<i>Peptide</i>	<i>%Cross-reactivity</i>
	Pro-BNP (1-45) (Rat)	100
	Pro-BNP-46 (Human)	0
	CNP-22 (Human, Porcine, Rat)	0
	BNP-45 (Rat)	0
	ANP (Rat)	0
Sequence	His - Pro - Leu - Gly - Ser - Pro - Ser - Gln - Ser - Pro - Glu - Gln - Ser - Thr - Met - Gln - Lys - Leu - Leu - Glu - Leu - Ile - Arg - Glu - Lys - Ser - Glu - Glu - Met - Ala - Gln - Arg - Gln - Leu - Ser - Lys - Asp - Gln - Gly - Pro - Thr - Lys - Glu - Leu - Leu	