

Catalog # RK-011-42

Standard Size 125 RIA tubes

Radioactivity This kit contains 1.5 µCi of Iodine 125

Range 10-1280 pg/tube in 100µl of Standard Solution = 100-12800 pg/ml

Lowest 602 pg/ml

Detection Limit

Cross Reactivity	Peptide	%Cross-reactivity
	N-terminal BNP (1-76) (Human)	100
	N-terminal BNP (47-76) (Human)	100
	Prepro -BNP (1-46) (Human)	0
	BNP 32 (Human)	0
	BNP32 (Rat)	0
	ANP-alpha (Human)	0
	CNP-22 (Human)	0
	Adrenomedullin (Human)	0
	Intermedin/AM-2 (Human)	0
	CGRP (Human)	0
	Angiotensin II (Human)	0
	Endothelin	0

Sequence HPLGSPGSASDLETSGLQEQRNHLQGKLSLQVEQTSLEPLQESPRTPGVWKS
REVATEGIRGHRKMVLYTLRAPR