

Catalog # RK-012-06

Standard Size 125 RIA tubes

Radioactivity This kit contains 1.5 µCi of Iodine 125

Range 10-1280 pg/tube in 100µl of Standard Solution =
100-12800 pg/ml

**Lowest
Detection Limit** 764.5 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	CNP-53 (1-29) (Porcine, Rat)	100
	CNP-53 (Porcine, Rat)	0
	CNP-22 (Human, Porcine, Rat)	0
	α-ANP (1-28) (Human, Canine)	0
	ANP (Rat)	0
	BNP-32 (Porcine)	0
	BNP-45 (Rat)	0
	Endothelin-1 (Human, Rat, Mouse, Porcine, Bovine, Canine)	0

Sequence Asp - Leu - Arg - Val - Asp - Thr - Lys - Ser - Arg - Ala - Ala
- Trp - Ala - Arg - Leu - Leu - His - Glu - His - Pro - Asn - Ala
- Arg - Lys - Tyr - Lys - Gly - Gly - Asn