

Catalog #	RK-015-02
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube of 100µl Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	27 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
-------------------------	----------------	--------------------------

CGRP (Human)	100
CGRP (8-37) (Human)	100
Amylin (Human)	0
Amylin Amide (Human)	0
ADM (Human)	0
Amylin (Rat)	0
Calcitonin (Salmon)	0
CRSP (24-38) Cys0 (Porcine)	0
CRSP (Porcine)	0
Substance P	0

Sequence Ala - Cys - Asp - Thr - Ala - Thr - Cys - Val - Thr - His - Arg - Leu - Ala - Gly - Leu - Leu - Ser - Arg - Ser - Gly - Gly - Val - Val - Lys - Asn - Asn - Phe - Val - Pro - Thr - Asn - Val - Gly - Ser - Lys - Ala - Phe - NH₂

Disulfide bridge between:Cys2 & Cys7