

Catalog #	RK-015-09
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube in 100µl of Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	47 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	CGRP (Rat)	100
	CGRP II (Rat)	78
	CGRP (Human)	35
	CGRP II (Human)	16
	Amylin (Rat)	0
	Amylin (Human)	0
	Calcitonin (Rat)	0
	Somatostatin-14	0

Sequence Ser - Cys - Asn - Thr - Ala - Thr - Cys - Val - Thr - His - Arg - Leu - Ala - Gly - Leu - Leu - Ser - Arg - Ser - Gly - Gly - Val - Val - Lys - Asp - Asn - Phe - Val - Pro - Thr - Asn - Val - Gly - Ser - Glu - Ala - Phe - NH₂

Disulfide bridge between: Cys2 & Cys7