

Catalog #	RK-017-01
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube in 100 µl of Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	50.49 pg/ml

Cross Reactivity	Peptide	%Cross reactivity
-------------------------	----------------	--------------------------

Amylin (Feline)	100
Amylin (Rat)	100
Amylin Amide (Human)	0
Amylin (Human)	0
CGRP (Human)	0
Insulin (Human)	0
Glucagon (Human)	0
Somatostatin-28	0
Pancreatic Polypeptide (Human)	0
C-Peptide (Human)	0

Sequence Lys - Cys - Asn - Thr - Ala - Thr - Cys - Ala - Thr - Gln - Arg - Leu - Ala - Asn - Phe - Leu - Ile - Arg - Ser - Ser - Asn - Asn - Leu - Gly - Ala - Ile - Leu - Ser - Pro - Thr - Asn - Val - Gly - Ser - Asn - Thr - Tyr - NH₂

Disulfide Bridge:
Between Cys2 & Cys7