

Catalog #	RK-019-26
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	10-2560 pg/tube in 100µl Standard Solution = 100-25600 pg/ml
Lowest Detection Limit	450 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Stresscopin (Human)	100
	Urocortin III (Human)	60
	Stresscopin related peptide (Human)	0
	Urocortin II (Human)	0
	Urocortin II (Mouse)	0
	Urocortin (Human)	0
	Urocortin (Rat)	0
	CRF (Human, Rat)	0

Sequence	Thr - Lys - Phe - Thr - Leu - Ser - Leu - Asp - Val - Pro - Thr - Asn - Ile - Met - Asn - Leu - Leu - Phe - Asn - Ile - Ala - Lys - Ala - Lys - Asn - Leu - Arg - Ala - Gln - Ala - Ala - Ala - Asn - Ala - His - Leu - Met - Ala - Gln - Ile - NH ₂
-----------------	---