

Catalog #	RK-021-03
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube in 100ul Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	23 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Dynorphin A (Porcine)	100
	Dynorphin A (1-13) (Porcine)	<1.0
	Dynorphin A (1-8) (Porcine)	0
	Dynorphin B (Porcine)	0
	β-Endorphin (Human)	0
	α-Neo-Endorphin (Porcine)	0
	Leu-Enkephalin	0

Sequence Tyr - Gly - Gly - Phe - Leu - Arg - Arg - Ile - Arg - Pro - Lys
- Leu - Lys - Trp - Asp - Asn - Gln