

Catalog #	RK-021-10
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube in 100µl Standard Solution = 10~1280 pg/ml
Lowest Detection Limit	23.12 pg/ml

Cross Reactivity	<i>Peptide</i>	<i>%Cross-reactivity</i>
	Dynorphin A (1-8) (Porcine)	100
	Big Dynorphin (Porcine)	0
	Dynorphin A (Porcine)	0
	Dynorphin A (1-13) (Porcine)	0
	Dynorphin B (Porcine)	0
	β-Endorphin (Porcine)	0
	α-Neo-Endorphin (Porcine)	0
	Leu-Enkephalin-Arg	0

Sequence Tyr - Gly - Gly - Phe - Leu - Arg - Arg - Ile