

Catalog #	RK-021-21
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube in 100µl Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	15.57 pg/ml

Cross Reactivity	<i>Peptide</i>	<i>%Cross-reactivity</i>
	Dynorphin A (1-13) (Porcine)	100
	Dynorphin A (Porcine)	100
	Dynorphin A (1-12) (Porcine)	100
	Big Dynorphin (209-240) (Porcine)	25
	Dynorphin A (1-10) Amide	1.5
	Dynorphin (1-11) (Porcine)	1.5
	Dynorphin (1-6, 7, 8, 9)	0
Sequence	Tyr - Gly - Gly - Phe - Leu - Arg - Arg - Ile - Arg - Pro - Lys - Leu - Lys	