

Catalog #	RK-022-14
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube of 100 µl Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	30.86 pg/ml

Cross Reactivity	<i>Peptide</i>	<i>%Cross-reactivity</i>
	β-Endorphin (Human)	100
	β-Endorphin (Rat)	0
	ACTH (Human)	0
	α-MSH	0

Sequence Tyr - Gly - Gly - Phe - Met - Thr - Ser - Glu - Lys - Ser - Gln - Thr - Pro - Leu - Val - Thr - Leu - Phe - Lys - Asn - Ala - Ile - Ile - Lys - Asn - Ala - Tyr - Lys - Lys - Gly - Glu