

Catalog #	RK-022-33
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube in 100 µl of Standard Solution = 10~1280 pg/ml
Lowest Detection Limit	27 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	β-Endorphin (Rat)	100
	β-Endorphin (Mouse, Porcine)	100
	β-Endorphin (Camel, Ovine, Bovine)	100
	β-Endorphin (Equine)	100
	α-Endorphin	100
	β-Endorphin (Human)	40
	Met-Enkephalin	0
	Leu-Enkephalin	0
	α-MSH	0
	ACTH (Human)	0
	PACAP 38 (Human, Ovine, Rat)	0

Sequence	Tyr - Gly - Gly - Phe - Met - Thr - Ser - Glu - Lys - Ser - Gln - Thr - Pro - Leu - Val - Thr - Leu - Phe - Lys - Asn - Ala - Ile - Ile - Lys - Asn - Val - His - Lys - Lys - Gly - Gln
-----------------	---