

Catalog #	RK-023-10
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube of 100µl Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	61.8 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Big Endothelin-1 (Human)	100
	Big Endothelin-1 (22-38) (Human)	100
	Big Endothelin-1 (Rat)	7
	Big Endothelin (22-39) (Porcine)	1
	Big Endothelin (Porcine)	0
	Endothelin-1 (Human, Rat, Porcine, Mouse, Canine, Bovine)	0
	Endothelin-2 (Human, Canine)	0
	Endothelin-3 (Human, Rat)	0
	BNP-32 (Porcine)	0
	Sarafotoxin S6b	0

Sequence Cys - Ser - Cys - Ser - Ser - Leu - Met - Asp - Lys - Glu - Cys - Val - Tyr - Phe - Cys - His - Leu - Asp - Ile - Ile - Trp - Val - Asn - Thr - Pro - Glu - His - Val - Val - Pro - Tyr - Gly - Leu - Gly - Ser - Pro - Arg - Ser

Disulfide Bridge: Cys1 - Cys15 and Cys3 - Cys11