

Catalog #	RK-023-13
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube of 100µl Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	603 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Endothelin-2 (Human)	100
	Endothelin-1 (Human, Rat, Porcine, Mouse, Canine, Bovine)	100
	Endothelin-3 (Rat)	100
	Big Endothelin-1 (Human)	100
	Big Endothelin-1 (Porcine)	100
	VIC (Mouse)	100
	Big Endothelin (22-38) (Human)	0
	Sarafotoxin S6b	0
	BNP 26 (Porcine)	0
	Angiotensin II (Human)	0
	Alpha-ANP (Human)	0

Sequence Cys - Ser - Cys - Ser - Ser - Trp - Leu - Asp - Lys - Glu - Cys - Val - Tyr - Phe - Cys - His - Leu - Asp - Ile - Ile - Trp

Disulfide Bridge: Cys1 - Cys15 and Cys3 - Cys11