

Catalog #	RK-023-37
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube in 100 µl of Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	60.44 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Prepro-Endothelin (18-50) (Human)	100
	Big Endothelin-1 (Human)	0
	Endothelin-1 (Human, Bovine, Canine, Mouse, Porcine, Rat)	0
	Endothelin-2 (Human, Canine)	0
	Endothelin-3 (Human, Porcine, Rat, Rabbit)	0

Sequence Ala - Pro - Glu - Thr - Ala - Val - Leu - Gly - Ala - Glu - Leu - Ser - Ala - Val - Gly - Glu - Asn - Gly - Gly - Glu - Lys - Pro - Thr - Pro - Ser - Pro - Pro - Trp - Arg - Leu - Arg - Arg - Ser