

Catalog #	RK-026-06
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube of 100µl Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	47 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Galanin (Porcine)	100
	Galanin (Human)	100
	Galanin (Rat)	100
	PHM-27 (Human)	0
	Secretin (Human)	0
	Insulin (Human)	0
	Substance P	0
	VIP (Human, Rat, Porcine)	0
	GMAP (1-41) Amide (Porcine)	0
	NPY (Porcine)	0

Sequence Gly - Trp - Thr - Leu - Asn - Ser - Ala - Gly - Tyr - Leu - Leu - Gly - Pro - His - Ala - Ile - Asp - Asn - His - Arg - Ser - Phe - His - Asp - Lys - Tyr - Gly - Leu - Ala - NH₂