

Catalog #	RK-026-52
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	10-1280 pg/tube of 100µl Standard Solution = 100-12800 pg/ml
Lowest Detection Limit	80.1 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	GALP (Rat)	100
	GALP (Human)	30
	GALP (36-60)	0
	Galantine (Rat)	0
	Galantine (Porcine)	0
	AGRP (82-131) (Mouse, Rat)	0
	AGRP (83-132) (Human)	0
	NPY	0
	α-MSH	0
	MCH	0
	CRF (Human)	0
	Somatostatin	0
	Orexin A	0
	Orexin B (Human)	0

Sequence	Ala - Pro - Ala - His - Arg - Gly - Arg - Gly - Gly - Trp - Thr - Leu - Asn - Ser - Ala - Gly - Tyr - Leu - Leu - Gly - Pro - Val - Leu - His - Leu - Ser - Ser - Lys - Ala - Asn - Gln - Gly - Arg - Lys - Thr - Asp - Ser - Ala - Leu - Glu - Ile - Leu - Asp - Leu - Trp - Lys - Ala - Ile - Asp - Gly - Leu - Pro - Tyr - Ser - Arg - Ser - Pro - Arg - Met - Thr -
-----------------	---