

Catalog # RK-027-02

Standard Size 125 RIA tubes

Radioactivity This kit contains 1.5 µCi of Iodine 125

Range 10-1280 pg/tube of 100µl Standard Solution =
100-12,800 pg/ml

**Lowest
Detection Limit** 376 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	GIP (Human)	100
	GIP (Porcine)	100
	GIP (Rat)	100
	Gastrin I (Human)	0
	Big Gastrin I (Human)	0
	Gastrin Releasing Peptide (Human)	0
	Bombesin	0
	Substance P	0
	Insulin (Human)	0
	Amylin Amide (Human)	0

Sequence Tyr - Ala - Glu - Gly - Thr - Phe - Ile - Ser - Asp - Tyr - Ser
- Ile - Ala - Met - Asp - Lys - Ile - His - Gln - Gln - Asp - Phe
- Val - Asn - Trp - Leu - Leu - Ala - Gln - Lys - Gly - Lys - Lys
- Asn - Asp - Trp - Lys - His - Asn - Ile - Thr - Gln