

Catalog # RK-028-14

Standard Size 125 RIA tubes

Radioactivity This kit contains 1.5 µCi of Iodine 125

Range 10-1280 pg/tube of 100µl Standard Solution =
100-12,800 pg/ml

**Lowest
Detection Limit** 1401 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	GLP-2 (1-34) (Human)	100
	GLP-2 (Rat)	100
	GLP-1 (Human, Rat, Mouse)	0
	GLP-1 (7-36) (Human, Rat, Mouse)	0
	GLP-1 (7-37) (Human, Rat, Mouse)	0
	Glucagon (Human)	0
	Insulin (Human)	0
	GIP (Human)	0
	VIP (Human, Porcine, Rat)	0

Sequence His - Ala - Asp - Gly - Ser - Phe - Ser - Asp - Glu - Met - Asn
- Thr - Ile - Leu - Asp - Asn - Leu - Ala - Ala - Arg - Asp -
Phe - Ile - Asn - Trp - Leu - Ile - Gln - Thr - Lys - Ile - Thr
- Asp - Arg