

Catalog # RK-028-39

Standard Size 125 RIA tubes

Radioactivity This kit contains 1.5 µCi of Iodine 125

Range 10-1280 pg/tube of 100µl Standard Solution =
100-12,800 pg/ml

**Lowest
Detection Limit** 153.3 pg/ml

Cross Reactivity	<i>Peptide</i>	<i>%Cross-reactivity</i>
	Pro-Resistin (43-65) (Human)	100

Sequence Ala - Ile - Ser - Ser - Ile - Gly - Leu - Glu - Cys - Gln - Ser
- Val - Thr - Ser - Arg - Gly - Asp - Leu - Ala - Thr - Cys -
Pro - Arg

Disulfide Bridge: Cys9 - Cys21