

Catalog #	RK-045-96
Standard Size	125 RIA Tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	4-512 pg/tube in 100µl of Standard Solution = 40-5120 pg/ml
Lowest Detection Limit	323 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Nm U, Prepro (104-136)(Human)	100
	NMS, Prepro (70-103)(Human)	0
	NMS, Prepro (70-103)(Rat)	0
	Neuromedin S (NMS) Amide (Human)	0
	Neuromedin S (Mouse)	0
	Nm U (Rat)	0
	Nm U-25 (Human)	0
	Ghrelin (Human)	0
	Motilin (Human)	0

Sequence	Phe - Leu - Phe - His - Tyr - Ser - Lys - Thr - Gln - Lys - Leu - Gly - Lys - Ser - Asn - Val - Val - Ser - Ser - Val - Val - His - Pro - Leu - Leu - Gln - Leu - Val - Pro - His - Leu - His - Glu
-----------------	---