

Catalog #	RK-065-32
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube of 100µl Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	35 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Copeptin (Human)	100
	Copeptin (Rat)	0
	Vasopressin (Human, Rat, Mouse, Canine)	0
	Angiotensin II (Human, Rat, Mouse, Canine)	0
	Angiotensin I (Human, Rat, Mouse, Porcine Canine)	0
	Endothelin-1 (Human, Rat, Mouse, Porcine, Bovine, Canine)	0
	Urotensin II (Human)	0

Sequence	Ala - Ser - Asp - Arg - Ser - Asn - Ala - Thr - Gln - Leu - Asp - Gly - Pro - Ala - Gly - Ala - Leu - Leu - Leu - Arg - Leu - Val - Gln - Leu - Ala - Gly - Ala - Pro - Glu - Pro - Phe - Glu - Pro - Ala - Gln - Pro - Asp - Ala - Tyr
-----------------	---